

SCOPE® Quattro

Ürün Özellikleri

SCOPE Quattro® 6 farklı kanal genişliği için 020-050 seçenekleri ile profesyonellerin kullanımına sunulmuştur. Sistemin, resiprokal hareketle kullanımı tavsiye edilir. İstek üzerine tek eğe sistemi olarak ta temin edilebilir. Daha esnek ve daha dayanıklı fiziksel özelliğiyle risksiz tedaviyi sağlar. Nitinol üzerindeki armor yapısı dayanımı artırırken mikron ölçüsü sayesinde, eğenin keskinlik ve esneklik performansını etkilememektedir.

*Eğeler Standart olarak 3 farklı boyutta 21,25,31mm olarak üretilebilmektedir.

SCOPE® Quattro

Product Feature

SCOPE Quattro is available to professionals with 6 different channel widths options of 020-050. It is recommended to use the system with reciprocal motion. Based on request, a single file system can be provided. It provides risk-free treatment with its more flexible and more durable physical features. The armor structure on nitinol increases the strength while the micron size does not affect the sharpness and flexibility performance of the file.

**The files can be produced in three different standard sizes of 21, 25, 31 mm.*



SCOPE[®] QUATTRO



Tipe Size / Uç çapı	Taper / Açı	Eğe Uzunluğu / File Length	Renk Kodu / Color Coding	Ürün Kodu / Product Code
20	07	21 / 25 / 31 mm	■	QG-(A)-(B)-(C)
25	07		■	
35	06		■	
45	05		■	
ASR1 (20-25-35-45)			■ ■ ■ ■	
(A)		(B)		Kesit / Cross Section
Paketleme / Packing	[1-6] BLISTER OR BOX (C)			



Gutta percha **GQG-(A)-(B)-(C)**
GQG-(20)-(06)-(25)
 Paper points **PQG-(A)-(B)-(C)**
PQG-(40)-(02)-(31)

Bu Modeller ile uyumlu çalışan gutta percha ve paper pointler
 Gutta percha and paper point for suitable with this models